

Gold Conjugate –Data Sheet

Conjugate:	Goat anti Rabbit IgG (H+L):	40nm Gold
Product Code:	BA.GAR40	
Lot Number:	21030154	

Description of test	Document reference	BBI Specification	Result	Pass/Fail
Optical Density at 520nm at manufacture:	QCP064	10.0 – 10.4	10.0	Pass
Mean Particle Diameter:	QCP007	37.0 – 43.0nm	40.5nm	Pass
Coefficient of Variation of Particle Size:	QCP007	≤ 8.0%	≤ 8.0%	Pass
Buffer:	Conjugate is suspended in 2mM sodium Tetraborate at pH 7.2 containing 0.095% Sodium Azide.			
Stability of this product is guaranteed when stored at 2-8°C until:	23 Mar 2023			

Ref: DS006 Issue 8
 Print Date: 30 Mar. 21
 Effective Date: Refer to QPulse